

# Innovation for people on the go

DreamStation Go provides the portability and packability travelers want and the effective therapy they need.



**Charge**  
Onboard USB port allows for easy charging of mobile devices.



**Power**  
Built-in power supply reduces the number of parts for you to manage for added convenience.



**Safe**  
Complies with FAA requirements for in-flight use.<sup>1</sup>



**Practical**  
Bluetooth is standard with every device.



**You're In control**  
Intuitive color touchscreen offers quick, simple setup and easy navigation.



**Compatibility**  
Device is compatible with our DreamMapper monitoring self-management tool.



<sup>1</sup>RTCA/DO-160G section 21, category M.



© 2017 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Caution: US federal law restricts these devices to sale by or on the order of a physician.

[www.philips.com/respironics](http://www.philips.com/respironics)

RRDPGH MK 3/14/17 MCI 4107896 PN 1134215  
1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 · 724-387-4000

**PHILIPS**

*DreamStation Go*



**Anywhere** you go

# DreamStation Go

Packs performance anywhere you go



## DreamStation Go

Whether you're traveling for business or venturing out on the vacation of a lifetime, DreamStation Go provides an easy, reliable and portable PAP experience when you refuse to compromise.

Our sleek design has 1/3 fewer components than other travel PAPs,<sup>1</sup> and at half the size of our trusted System One and DreamStation platforms, DreamStation Go offers the same therapy and comfort features enjoyed by over 5 million PAP users.

<sup>1</sup>When compared to against HDM Z1 and Somnetics Transcend devices (2016)



## Convenience without compromise

Wherever your travels take you, sleep apnea follows. That's why DreamStation Go is an indispensable travel partner, packing powerful performance into a light, streamlined package. It even has a convenient USB port to charge your mobile device.

- Our smallest and lightest tubing<sup>1</sup>
- Our smallest PAP machine
- Our lightest PAP at less than two pounds
- Complies with FAA requirements<sup>2</sup>

<sup>1</sup>When compared to Philips 22mm performance tubing  
<sup>2</sup>RTCA/DO-160G section 21, category M.

## CPAP therapy even 'off the grid'

### Overnight battery

When traveling to locations off the power grid, it can be difficult for those with sleep apnea to maintain their therapy. That's why we developed the DreamStation Go overnight battery.



DreamStation Go set at 10cm of pressure displayed an average run time of 13 hours<sup>1</sup> when running off the DreamStation Go overnight battery.

<sup>1</sup>PAP run-time test conditions: CPAP mode, pressure 10 cm H<sub>2</sub>O, 12mm tubing, 37LPM of leak, 73.4 F room at an elevation of 1650 ft.



### Small travel kit

Kit includes an easy-to-pack protective case for the DreamStation Go and a separate, washable bag for your mask and tube.



### Medium travel kit

Kit includes a protective case that accommodates the DreamStation Go and overnight battery. Also includes separate, washable bag for your mask and tube.

## Learn more

Visit [sleepapnea.com/dreamstationgo](http://sleepapnea.com/dreamstationgo) or call 1-844-240-1649.